HipFix® Thermoplastic Positioning System

CIVCO RT Technical Data Sheet



Product Name	HipFix Thermoplastic Positioning System. Includes baseplate and thermoplastic options.			
Product Number	MTPLVC04MR HipFix Baseplate, indexable, 2 & 3-Pin Lok-Bar compatible See page 3 for thermoplastics			
General Description	ne HipFix hip and pelvis immobilization system utilizes a single sheet of thermoplastic over e entire abdomen and pelvis region, providing immobilization to enhance fixed fiducial rgeting. The versatile system accommodates both prone and supine patient setups.			
Indications for use	Device is intended for patient fixation. Devices are used to aid in the support and positioning of patients during an MRI.			
MRI Compatible	MR Safe			
Integration	2 or 3-Pin Indexable to couchtop			
FDA Clearance	Yes			
CE Marked (EU DOC)	Yes			
Other Country Registrations	Available upon request			
Country of Manufacture	USA			




Baseplate Technical Details

Primary Materials	Black PVC, Gray PVC, White Delrin				
Material Density (approx.)	1.42 g/cm ³				
Weight (approx)	4.6 kg / 10.15 lbs				
Dimensions (mm)	787 L x 603 W x 49 H (min.)				
	Location*	Field Size (cm²)	Attenuation 6 MV photons	Attenuation 18 MV photons	
Dosimetric Properties	Please note: Imaging distortions an	6 x 6 d system reported Houns	5.19% field Units variations are	3.61% possible at the superior or	

Please note: Imaging distortions and system reported Hounsfield Units variations are possible at the superior or inferior ends of the HipFix based on the CT scanner or imaging settings. While the board is manufactured from the same material in these regions, it is not recommended to treat through these areas without making proper adjustments in the imaging or planning system.

Performance Specifications are for reference only, values may vary based on setup and equipment.

Technical Drawing (MTPLVC04MR)





Thermoplastic Technical Details

IMRT Reinforced™					
	Thickness	Weight	Dimensions	0 0 0	
MTHFI1822 (reloadable)	3.2 mm	1 lb (0.45 kg)	559 L x 457 W		
MTHFI18222.4 (reloadable)	2.4 mm	1 lb (0.45 kg)	559 L x 457 W		
MTPLHAI0203 (no groin piece, disposable)	3.2 mm	1.2 lb (0.54 kg)	559 L x 445 W	Reloadable	l : : : : : : : : : : : : : : : : : : :

Perforated			
	Thickness	Weight	Dimensions
MTHFX1822 (reloadable)	3.2 mm	1 lb (0.45 kg)	559 L x 457 W
MTHFX18222.4 (reloadable)	2.4 mm	1 lb (0.45 kg)	559 L x 457 W
MTPLHA0202 (no groin piece, disposable)	3.2 mm	1.76 lb (0.80 kg)	559 L x 461 W

Solid					
	Thickness	Weight	Dimensions	0 0 0	
MTHFX1822S (reloadable)	3.2 mm	1 lb (0.45 kg)	559 L x 457 W		4
MTHFX18222.4S (reloadable)	2.4 mm	1 lb (0.45 kg)	559 L x 457 W		۰
MTPLHA0201 (no groin piece, disposable)	3.2 mm	1.2 lb (0.54 kg)	559 L x 445 W	Reloadable	Disposabl



Thermoplastic Technical Details

2.4 Standard Perforation 3.2 Standard Perforation 2.4 IMRT Reinforced 3.2 IMRT Reinforced

Primary Materials	Polycaprolactone, White ABS							
Attenuation (Unstretched Mask)**	6MV: 0.53% 18MV: 0.03		6MV: 0.72% 6MV: 0.81% (solid region) 18MV: 0.12% 18MV: 0.21% (solid region)			6MV: 1.09% (solid region) 18MV: 0.21% (solid region)		
Thickness	2.4 mm (.09	94")	3.2 mm (.12	25")	2.4 mm (.09	4")	3.2 mm (.12	25")
Perforation (surface)	38% (Stand	ard)	38% (Stand	ard)	Varies (IMR	T Pattern)	Varies (IMRT Pattern)	
Density				1.145	g/cm3			
Surface feeling				Sm	ooth			
Color				W	nite			
Odor	None							
Biocompatible	Yes							
Transparent when activated				Y	es			
	Waterbath	Dry Heat Oven	Waterbath	Dry Heat Oven	Waterbath	Dry Heat Oven	Waterbath	Dry Heat Oven
Activation temperature	75° C / 165° F							
Activation time	6 min 4–8 min		6 min	4–8 min	7–9 min	8–11 min	7–9 min	8–11 min
Molding time (maximum)*a	1 min 15 sec	N/A	1 min 30 sec	N/A	1 min 30 sec	N/A	1 min 30 sec	N/A
Hardening time*b	8 min 45 sec	N/A	8 min 30 sec	N/A	8 min 30 sec	N/A	8 min 30 sec	N/A
Time to completion	10 min 10 min		10 min	10 min	10 min	10 min	10 min	10 min

Definitions*

- a) Molding Time: the time frame available from when the mask is removed from the waterbath to the time that the mask can no longer be molded or formed.
- b) Hardening Time: the time between the molding and completion time where the mask should be cooling and hardening on its own without being manipulated or molded.
- c) **Time to Completion**: the total time from when the mask is taken out of the waterbath or dry heat oven to when the mask can be removed from the patient.

Attenuation Test Conditions**

- a) 100 cm SAD
- b) 10 x 10 cm Field Size
- c) 6 MV Dmax = 1.5 cm
- d) 18 MV Dmax = 3.0 cm



Packaging

Unit Pack	Each		
Units / Box	Baseplate: One device – may have several component parts in box Thermoplastic: Qty 5		
Shelf Life	n/a		
Warranty	See CIVCO RT website for Product Warranty Statement		
Sterilization	Non-Sterile		
Launch Date	July 2019		

Accessories

	Baseplate Compatible with:				
System Customization	MTHFI1822 MTHFI18222.4 MTHFX1822 MTHFX18222.4 MTHFX1822S MTHFX18222.4S MTPLHA0201	IMRT Reinforced Thermoplastic, reloadable, 3.2 mm IMRT Reinforced Thermoplastic, reloadable, 2.4 mm Perforated Thermoplastic, reloadable, 3.2 mm Perforated Thermoplastic, reloadable, 2.4 mm Solid Thermoplastic, reloadable, 3.2 mm Solid Thermoplastic, reloadable, 2.4 mm Solid Thermoplastic, reloadable, 2.4 mm Solid Thermoplastic, disposable (no groin piece), 3.2 mm			
	MTPLHA0202 MTPLHAI0203 MTPLVCSUB8 MTPLHA01	Perforated Thermoplastic, disposable (no groin piece), 3.2 mm IMRT Reinforced Thermoplastic, disposable (no groin piece), 3.2 mm HipFix disposable thermoplastic groin piece HipFix Reloadable Handles (x2) and Groin Piece			

