Body Pro-Lok ONE™ Rails-Only for Varian Exact® and IGRT Style Couchtops

CIVCO RT Technical Data Sheet




Product Name	Body Pro-Lok ONE Varian Exact/IGRT Rails- Only
Product Number	BPL020
Kits that include this product	BPL300
General Description	Allows for connection of Body Pro-Lok ONEBridge™, ONE Respiratory Belt and other Body Pro-Lok ONE devices to aid in providing comfortable, total body immobilization for hypo-fractionated treatments
Indications for use	SBRT Patient Immobilization
MRI Compatible	No
Integration	Clamps to the IGRT couchtop profile and Exact Indexing CT overlays with an edge thickness of 0.76"-0.975" (19.3 mm-24.8 mm)
FDA Clearance	Yes
CE Marked (EU DOC)	Yes
Other Country Registrations	Available upon request



Technical Details

Primary Materials	Carbon Fiber, Black Acetal, Aluminum, Stainless Steel
Weight (approx)	2.49 kg / 5.5 lbs
Dimensions (mm)	367 L x 617 W x 81 H
Attenuation	Carbon Fiber Crossbar: 6 MV: 4.6% 18 MV: 3.0%

Technical Drawing



ry of Manufacture USA	y of Manufacture USA	Country of Manufacture US
-----------------------	----------------------	---------------------------



Packaging

Unit Pack	Each
Units / Box	One device - may have several component parts in box
Shelf Life	n/a
Warranty	See CIVCO RT website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	November 2016

Compatible Accessories

System Customization ONEBridge™: BPL001

ONE Respiratory Belt: BPL026

