


Bolus

CIVCO Technical Data Sheet



Product Name	Bolus with Skin, Skinless Bolus		
Product Number	Skinless (Not CE Marked)	Skinless (CE Marked)	With Skin (CE Marked)
MTCB405		MTCB403	MTCB403S
MTCB410		MTCB403CE	MTCB405S
MTCB415		MTCB405CE	MTCB410S
MTCB420		MTCB410CE	MTCB415S
		-	MTCB420S
-	-	-	MTCB305005S
-	-	-	MTCB305010S
-	-	-	MTCB450S
-	-	-	MTCB455S
-	-	-	MTCB404005S
-	-	-	MTCB404010S

General Description

Bolus is made from a solid, homogeneous, tissue-equivalent gel. It will not sag, deform or run. Bolus is available encased in a tough layer of thin plastic “skin” or “skinless”, which is sticky for greater conformity to steep body contours. Bolus sheets can be cut with a scissors to any desired shape and layered to obtain the desired bolus effect.

Indications for use	Bolus is intended as a tissue equivalent gel pad.
CE Marked (EU DOC)	Bolus w/ Skin: Yes Skinless Bolus: Yes (MTCB403, MTCB403CE, MTCB405CE, MTCB410CE)
Other Country Registrations	Available upon request
Country of Manufacture	USA

Technical Details

Primary Materials	Solid, homogenous, tissue-equivalent gel. Skin: Polyurethane/Synthetic Rubber		
Weight (approx) + Dimensions	Part Number	Weight (gm)	Dimensions
	MTCB403S	278.1	30 x 30 x 0.3 cm
	MTCB403	275.4	30 x 30 x 0.3 cm
	MTCB403CE	278.1	30 x 30 x 0.3 cm
	MTCB405S	463.5	30 x 30 x 0.5 cm
	MTCB405	748.4	30 x 30 x 0.5 cm
	MTCB405CE	459	30 x 30 x 0.5 cm
	MTCB410S	927	30 x 30 x 1.0 cm
	MTCB410	918	30 x 30 x 1.0 cm
	MTCB410CE	918	30 x 30 x 1.0 cm
	MTCB415S	1390.5	30 x 30 x 1.5 cm
	MTCB415	2177.2	30 x 30 x 1.5 cm
	MTCB420S	1854	30 x 30 x 2.0 cm
	MTCB420	2812.3	30 x 30 x 2.0 cm
	MTCB305005S	772.5	30 x 50 x 0.5 cm
	MTCB305010S	1545	30 x 50 x 1.0 cm
	MTCB455S	115.9	15 x 15 x 0.5 cm
	MTCB450S	231.8	15 x 15 x 1.0 cm
	MTCB404005S	824	40 x 40 x 0.5 cm
	MTCB404010S	1648	40 x 40 x 1.0 cm

Packaging

Unit Pack	1 piece
Units / Box	1 piece
Shelf Life	n/a
Warranty	See CIVCO website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	May 2006

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

