

# C-Qual M™ Breastboard with Monarch™ Overhead Arm Positioner

## CIVCO Technical Data Sheet



<b>Product Name</b>	C-Qual M™ Breastboard	
<b>Product Number</b>	MTM400	C-Qual M Supine Breastboard with Monarch Overhead Arm Positioner
	OAP100	Monarch Overhead Arm Positioner for use with C-Qual™
	OAP110	Monarch Overhead Arm Positioner for use with Lok-Bar™
	OAP120	Monarch Thermoplastic Frame Only, Upgrade Kit for Monarch Overhead Arm Positioner
<b>General Description</b>	The C-Qual M assists in providing comfortable, reproducible patient positioning for supine breast treatment. The board features an enhanced bottom stop. The board integrates with the Monarch Overhead Arm Positioner to provide an efficient approach for arms-up positioning, allowing for proper superior/inferior adjustment of the board to the patient which helps reduce patient slide.	
<b>Indications for use</b>	OAP110	The device is intended to support and aid in positioning a patient during radiologic and other medical procedures. Devices are used to aid in the support and positioning of patients during an MRI.
	MTM400 including OAP100	The device is indicated to aid in supporting and positioning adult and adolescent patients undergoing radiation therapy of the breast and chest region including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning. The device is not intended for use with patients under 12 years of age.
<b>MRI Compatible</b>	OAP110	MR Conditional. Product poses no known hazard in 1.5 and 3.0T MR environments.
<b>Integration</b>	3-Pin Indexable to Couchtop	
<b>FDA Clearance</b>	Yes	
<b>CE Marked (EU DOC)</b>	Yes	
<b>Other Country Registrations</b>	Available upon request	

# Technical Details

**Primary Materials** Kevlar, Black Acetal, Grey Kydex, Black ABS


**Weight (approx)** 10.89 kg / 24 lbs

**Dimensions (mm)** 1333 L x 594 W x 261 H (min.)

Dosimetric Properties	Location	6 MV photons	Proton water equivalency R80
	Middle of the back support ●	1.87	–
	Right Side ●	–	3.5 mm

Performance Specifications are for reference only, values may vary based on setup and equipment.

## Technical Drawing



**Country of Manufacture** USA

---

# Packaging

---

<b>Unit Pack</b>	Each
<b>Units / Box</b>	One device - may have several component parts in box
<b>Shelf Life</b>	n/a
<b>Warranty</b>	<a href="#">See CIVCO website for Product Warranty Statement</a>
<b>Sterilization</b>	Non-Sterile
<b>Launch Date</b>	October 2015

---

# Accessories

---

<b>System Customization</b>	Compatible with: Posicast® Head-Only Thermoplastics Mastercast™ Pro Head-Only Thermoplastics Posifix® Headrests Breast Support Thermoplastics Vac-Lok™ Cushions C-Qual Storage: 013156
-----------------------------	--

---

## Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

