

# Chabner XRT® Radiation Bra

## CIVCO Technical Data Sheet



<b>Product Name</b>	Chabner XRT® Radiation Bra			
<b>Product Number</b>	<b>CRTBRA01</b>	Size 1, 34"–36", B/C		
	<b>CRTBRA02</b>	Size 2, 34"–36", D/DD		
	<b>CRTBRA03</b>	Size 3, 38"–40", B/C		
	<b>CRTBRA04</b>	Size 4, 38"–40", D/DD		
	<b>CRTBRA05</b>	Size 5, 38"–40", F/G		
	<b>CRTBRA06</b>	Size 6, 42"–44", D/DD		
	<b>CRTBRA07</b>	Size 7, 42"–44", F/G		
	<b>CRTBRA08</b>	Size 8, 46"–48", D/DD		
	<b>CRTBRA09</b>	Size 9, 46"–48", F/G		
	<b>CRTEXT01</b>	Chabner XRT Radiation Bra Extender (one size)		
<b>Kits that include this product</b>	CRTBRASPK	CRTBRA01PK5	CRTBRA04PK5	CRTBRA07PK5
	CRTBRASPKLT	CRTBRA02PK5	CRTBRA05PK5	CRTBRA08PK5
	CRTEXT01PK5	CRTBRA03PK5	CRTBRA06PK5	CRTBRA09PK5
<b>General Description</b>	Provides optimal breast support, patient comfort and dignity, and reproducibility during simulation, planning and radiation treatment			
<b>Indications for use</b>	This device is intended to aid in the support of the patient during X-ray, computed tomography, magnetic resonance imaging, radiotherapy, and other diagnostic radiological procedures. Not intended for pediatric use.			
<b>MRI Compatible</b>	Yes			
<b>Integration</b>	n/a			
<b>FDA Clearance</b>	Yes			
<b>CE Marked (EU DOC)</b>	No			
<b>Other Country Registrations</b>	n/a			
<b>Country of Manufacture</b>	USA			

---

# Technical Details

---

## Primary Materials

UBL Tricot, Mesh Warp Knit, Synthetic Elastomer, Nylon, Polyester

---


## Attenuation

Material	6 MV photons	6 MeV electrons
Tricot	0%	0%
TPU	0%	0%
Elastic	0%	.1%
Velcro	0%	.7%
Mesh	0%	0%

Performance Specifications are for reference only, values may vary based on setup and equipment.

---

## Technical Drawing



Part Number	Cup Size	Dimensions: Total Bottom Band on the Half (Fully Closed)	Product Weight (Approx)
CRTBRA01	B/C	14.25 in (36.2 cm)	2.1 oz (60 grams)
CRTBRA02	D/DD	14.25 in (36.2 cm)	2.3 oz (65 grams)
CRTBRA03	B/C	15.75 in (40 cm)	2.5 oz (70 grams)
CRTBRA04	D/DD	15.75 in (40 cm)	2.6 oz (75 grams)
CRTBRA05	F/G	16.00 in (40.6 cm)	2.8 oz (80 grams)
CRTBRA06	D/DD	17.5 in (44.5 cm)	3.0 oz (85 grams)
CRTBRA07	F/G	17.375 in (44.1 cm)	3.2 oz (90 grams)
CRTBRA08	D/DD	18.875 in (47.9 cm)	3.4 oz (95 grams)
CRTBRA09	F/G	18.75 in (47.6 cm)	3.5 oz (100 grams)

## Packaging

<b>Unit Pack</b>	Each
<b>Units / Box</b>	One device
<b>Shelf Life</b>	n/a
<b>Warranty</b>	<a href="#">See CIVCO website for Product Warranty Statement</a>
<b>Sterilization</b>	Non-Sterile
<b>Launch Date</b>	October 2017

### Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

