

CIVCO Technical Data Sheet

MT-350 Breastboard



Product Name	MT-350 Breastboard
Product Number	MT350N
General Description	The MT-350 assists in providing comfortable, reproducible patient positioning for supine breast treatment. The breast board features a positive-locking tilt system, helping to ensure proper breast set-up angulation. Its carbon fiber treatment grid and adjustable arm and wrist supports aid in comfortable arms-up positioning for enhanced breast treatment access.
Indications for use	Breast cancer patients unable to lie prone. Smaller-breasted patients. Need for clearer visualization of medial field edge. Enhanced CT-gantry clearance.
MRI Compatible	No
Integration	2-Pin Indexable to Couchtop
FDA Clearance	Yes
CE Marked (EU DOC)	Yes
Other Country Registrations	Available upon request


Technical Details

Primary Materials Carbon Fiber

Weight (approx) 6.62 kg / 14.6 lbs

Dimensions (mm) 1257 L x 463 W x 92 H (min.)

Technical Drawing



Country of Manufacture USA

Packaging

Unit Pack	Each
Units / Box	One device - may have several component parts in box
Shelf Life	n/a
Warranty	See CIVCO website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	June 1998

Accessories

System Customization	Round Head Cup: MT250SUB4 Gel Head Cup: MT250SUB5 Timo Headrest Adapter Plate: MT350400 Right Arm Support: MTTASR Left Arm Support: MTTASL Wrist Support: MTTWSS Delrin Hand Grips: MT2504 MT-350 Wall Storage Device: MT350WSD
-----------------------------	--

Committed to Positively Impacting Patient Care

Global Offices
Orange City, Iowa, USA Rotterdam, Netherlands info@zivco.com
800.842.8688 | +1.712.737.8688 +31 (0) 10.3036600 www.zivco.com/ro

COPYRIGHT © 2015, CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2015C1011 REV B

