


Universal Couchtop™ 1 Piece

CIVCO Technical Data Sheet




| | | |
|---------------------------------------|--|--|
| Product Name | Universal Couchtop™ 1 Piece | |
| Product Number | CT: Siemens Biograph (Short): MTIL6022 Siemens Somatom (Long): MTIL6023 Siemens Definition (Long): PHS 1 A/B Cradle MTIL6024 PHS 4 Cradle MTIL6026 GE HiSpeed (Short): MTIL6320 GE Discovery ST (Long): MTIL6321 GE Global 1700 (Short): MTIL6322 GE Global 2000 (Long): MTIL6323 Philips Brilliance (Long): MTIL6420 Toshiba Aquilion (Short): MTIL6520 Toshiba Aquilion One (Long): MTIL6521 | Linac: Siemens Txt (Long): MTIL6020 Siemens ZXT (Long): MTIL6021 Varian Exact (Long): MTIL6120 Protura on Varian (Long): MTIL6121 Elekta (Long): MTIL6220 |
| Kits that include this product | Not Applicable | |
| General Description | The Universal Couchtop 1 Piece is a strong and lightweight couchtop with low attenuation and surface dose. It is available for most CT scanners and linear accelerators, providing a dosimetrically matched treatment platform. | |
| Indications for use | To support and aid in positioning a patient during radiologic and other medical procedures. | |
| Contraindications | None | |
| MRI Compatible | No | |
| Integration | Not Applicable | |
| FDA Clearance | Yes | |
| CE Marked (EU DOC) | Yes | |
| Other Country Registrations | Available upon request | |

Technical Details


| | Long | Short |
|-------------------------------|------------------------|------------------------|
| Primary Materials | Carbon Fiber Composite | Carbon Fiber Composite |
| Weight (approx) | 13.6 kg / 30 lbs | 12.3 kg / 27 lbs |
| Dimensions (mm) | 2432 L x 530 W x 50 H | 2165 L x 530 W x 50 H |
| Country of Manufacture | USA | USA |

Technical Drawing

Long



Short



Packaging

| | |
|----------------------|--|
| Unit Pack | 1 Piece |
| Units / Box | 1 Piece |
| Shelf Life | Not Applicable |
| Warranty | See CIVCO website for Product Warranty Statement |
| Sterilization | Non-Sterile |
| Launch Date | June 2010 |

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

