


# Universal Couchtop™ 1 Piece

## CIVCO Technical Data Sheet




<b>Product Name</b>	Universal Couchtop™ 1 Piece	
<b>Product Number</b>	<b>CT:</b> Siemens Biograph (Short): MTIL6022 Siemens Somatom (Long): MTIL6023 GE HiSpeed (Short): MTIL6320 GE Discovery ST (Long): MTIL6321 GE Global 1700 (Short): MTIL6322 GE Global 2000 (Long): MTIL6323 Philips Brilliance (Long): MTIL6420 Toshiba Aquilion (Short): MTIL6520 Toshiba Aquilion One (Long): MTIL6521	<b>Linac:</b> Siemens Txt (Long): MTIL6020 Siemens ZXT (Long): MTIL6021 Varian Exact (Long): MTIL6120 Protura on Varian (Long): MTIL6121 Elekta (Long): MTIL6220
<b>Kits that include this product</b>	Not Applicable	
<b>General Description</b>	The Universal Couchtop 1 Piece is a strong and lightweight couchtop with low attenuation and surface dose. It is available for most CT scanners and linear accelerators, providing a dosimetrically matched treatment platform.	
<b>Indications for use</b>	To support and aid in positioning a patient during radiologic and other medical procedures.	
<b>Contraindications</b>	None	
<b>MRI Compatible</b>	No	
<b>Integration</b>	Not Applicable	
<b>FDA Clearance</b>	Yes	
<b>CE Marked (EU DOC)</b>	Yes	
<b>Other Country Registrations</b>	Available upon request	

# Technical Details


	Long	Short
<b>Primary Materials</b>	Carbon Fiber Composite	Carbon Fiber Composite
<b>Weight (approx)</b>	13.6 kg / 30 lbs	12.3 kg / 27 lbs
<b>Dimensions (mm)</b>	2432 L x 530 W x 50 H	2165 L x 530 W x 50 H
<b>Country of Manufacture</b>	USA	USA

## Technical Drawing

Long



Short



---

# Packaging

---

<b>Unit Pack</b>	1 Piece
<b>Units / Box</b>	1 Piece
<b>Shelf Life</b>	Not Applicable
<b>Warranty</b>	<a href="#">See CIVCO website for Product Warranty Statement</a>
<b>Sterilization</b>	Non-Sterile
<b>Launch Date</b>	June 2010

---

Global Sales Office  
Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@civco.com  
www.civco.com

COPYRIGHT © 2017, CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. UNIVERSAL COUCHTOP AND PROTURA ARE TRADEMARKS OF CIVCO.  
ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.  
2015C0966 REV. D

**CIVCO**<sup>®</sup>  
Radiotherapy