

# Hip + Pelvic Positioning Devices


Radiation therapy treatments in the pelvic region often utilize fields with very small margins. This makes accurate patient positioning an essential step in the treatment process. An accurate, consistently reproducible patient setup throughout the entire course of treatment can only be achieved with the right equipment. CIVCO offers a variety of products which can be used to treat patients in either the prone or supine position.

## Carbon Fiber Bellyboard

The Carbon Fiber Bellyboard was developed to reduce the volume of the irradiated small bowel for patients undergoing treatment in the pelvic region.

- Fiducial markers integrate in the bottom part of the device provide **patient setup verification**
- Clear scales on both sides of the bellyboard allow **easy verification of the longitudinal setup**
- Radiotranslucent making it **radiation-friendly**
- Lightweight and rigid providing clinical **ease of use**
- Indexable for **quick, accurate setup**

A Prone Thorax Support is available to relieve chest compression during treatments in the prone position. It can be used as a stand alone patient support or as an option on the bellyboard to allow large patients to be treated in the prone position.



Carbon Fiber Bellyboard	125012
MR Bellyboard	126000
Prone Thorax Support	301102
<i>sold separately</i>	

## Combifix™

The Combifix is a baseplate which combines the Kneefix™ 3 and Feetfix™ 3. It is recommended both cushions are used together, as each cushion enhances the performance of the other, especially when treating the pelvic region. The Combifix provides a wide variety of adjustment possibilities, including a shoulder retractor option, providing a **custom setup for each patient**.

The Kneefix 3 cushion contributes to better stability and reproducibility in the supine position, regardless of the target area. It promotes patient comfort and is therefore ideal for use during lengthy procedures on rigid surfaces.

The Feetfix 3 is a universal patient support cushion for radiotherapy and diagnostic procedures and has been designed to provide **comfortable positioning with enhanced immobilization**.



Combifix (includes 1 Combifix baseplate, 1 Kneefix, 2 Kneefix Elevation Blocks, 1 Feetfix and 1 Feetfix Couch Support)	151002
Kneefix	301052
Kneefix Elevation Block	301055
Kneefix Elevation Block (set of 3)	301056
Feetfix	301071
Shoulder Retractors	302070

Contact your CIVCO representative for more information

[WWW.CIVCO.COM/RO](http://WWW.CIVCO.COM/RO) | 800.842.8688 | +1 712.737.8688

COPYRIGHT © 2015, CIVCO IS A REGISTERED TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. COMBIFIX, KNEEFIX, AND FEETFIX ARE TRADEMARKS OF CIVCO. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2015P0973 REV.A



Radiation Oncology