

Universal Couchtop™ for Safe, Efficient Treatment Planning & Delivery

 ProtonSeries™

Junction, Rail and Artifact-free

Specifically designed to answer the demands of Proton, the Universal Couchtop Long Extension contains no rails or junctions in the treatment area to help facilitate better dose delivery and target coverage.

Universal Couchtop™ Long Extension

- Large treatment area of 64.9" (164.8 cm)
- Homogenous design for consistent treatment planning & delivery
- Utilizes Prodigy™ 2 indexing at 7 cm increments
- Low attenuation and surface dose
- Self-locking extension seam located outside of treatment field
- 5 Year Warranty



IBA Systems Proteus® ONE and Proteus PLUS

MTIL6840 F-Series

MTIL6841 L-Series

Both part numbers include: UCT LE Couchtop with Rectangular Extension (MTIL6673) and Prodigy™ 2 Lok-Bar™



Global Sales Office
Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com



CIVCO®
Radiotherapy

COPYRIGHT © 2018, CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. PRODIGY AND UNIVERSAL COUCHTOP ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2018P1354 REV. B