## Confident Pediatric Proton Therapy without Compromise



**ProForm**™

The Universal Couchtop™ provides confidence in your pediatric care with no rails or junctions in the treatment area to facilitate better access, delivery of dose and target coverage.

- Large, 108.6 cm artifact-free treatment area with quick and easy attachment
- Extension features geometry that allows for direct access to effectively treat CNS pediatric patients
- Homogenous treatment area and soft gradient edges reduce beam degradation
- Prodigy<sup>™</sup> 2 indexing with 7 cm spacing
- Thermoplastic mask attaches with a smaller profile frame and pin similar to CIVCO RT's Type-S™ system





NEW

MTIL6682K

Pediatric ProForm Head & Neck Extension, 48.7" (123.7 cm), includes Cushion Guide and UCT Patient Handles



## Pediatric ProForm™ Head & Neck Overlays

NEW

The ProForm CT & Photon Treatment Overlay and MR Overlay allow for CT simulation or MR imaging of patients who will be treated on the Universal Couchtop ProForm Extension, allowing seamless transition from simulation to treatment.

PRO400K Pediatric ProForm CT & Photon Treatment Overlay,

2-Pin Lok-Bar compatible, includes ProForm Cushion

Guide and Patient Handles

Pediatric ProForm MR Overlay, 3-Pin Lok-Bar compatible, PRO420K

includes ProForm Cushion Guide and Patient Handles



## Thermoplastic Styles





PRO430IM1832CL Pediatric Head Only

PRO431IM1832CL Pediatric Head. Neck & Shoulders

ClearVision™ 2, Disposable, 3.2 mm

PRO430IM1832SG Pediatric Head Only

PRO431IM1832SG Pediatric Head. Neck & Shoulders



PRO430IM1832 Pediatric Head Only

PRO431IM1832 Pediatric Head, Neck & Shoulders

Standard Perforated Thermoplastic, Disposable, 3.2 mm

PRO430 Pediatric Head Only















## Visit CivcoRT.com for more information!









info@CivcoRT.com | www.CivcoRT.com

