

Corrective Pitch Capability & Positioning Flexibility



Solstice™
Immobilization System

Solstice™ Immobilization System

CIVCO Radiotherapy brings leading-edge innovation and advancement of non-invasive, head and neck treatment with the Solstice Immobilization System. The system provides unique corrective pitch capability and positioning flexibility.

The system is comprised of a carbon fiber head support, customizable cushion and dedicated thermoplastic mask designed for secure and simple attachment. These components create a highly customizable and comfortable head and neck immobilization system for precise treatment to improve patient outcomes.

510(k) pending



Featuring IMRT Reinforced™ Thermoplastic



Featuring ClearVision™ Open-face Thermoplastic with Precise Bite™, compatible with VisionRT's AlignRT® camera system

Features & Clinical Benefits

- Variable pitch capability provides increased setup options in CT simulation and corrective positioning flexibility during treatment setup
- Midline thermoplastic mask attachment provides increased rigidity and minimizes shrinkage to enhance patient comfort
- Highly conformal to the patient when using a customizable cushion and dedicated thermoplastic mask
- Adaptable with Precise Bite to aid in effective and repeatable positioning
- Unique thermoplastic frame provides easy to use, secure and quiet mask attachment
- Type-S™ and Posifix® version compatible



Pitch at 0 (chin down)



Pitch at 5 (neutral)



Pitch at 10 (chin up)

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

COPYRIGHT © 2018. CIVCO AND POSIFIX ARE REGISTERED TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. CLEARVISION, IMRT REINFORCED, PRECISE BITE, SOLSTICE AND TYPE-S ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2018P1357 REV. B

CIVCO®
Radiotherapy