

**NEW!**

# Vac-Lok™ Solutions for Monarch™ Overhead Arm Positioner & Wing Boards™

CIVCO offers a full-line of Vac-Lok cushions for head & neck, breast, hip & pelvic and full body positioning needs. Vac-Lok Cushions create a rigid and secure support around the patient when a vacuum is drawn through a self-sealing quick-release valve.

We are pleased to share two new Vac-Lok cushion offerings that provide additional size and molding fill options for positioning setups when used with the Monarch Overhead Arm Positioner or other CIVCO Wing Boards.

The new cushions retain the current features and benefits of all Vac-Lok, including:

- Retain shape and guarantee stability and reproducibility of the patient's position for up to six weeks once the air is vacuumed out
- Protective material coating provides greater cleanability than previous generation bags
- Treatment area is artifact-free with minimal beam attenuation
- Reusable for ultimate cost-effectiveness

In addition to these points, the new cushions feature these design enhancements:

- Offers a larger cushion to capture a greater amount of the patient's arms (compared VLWGB03, 88 x 55 cm)
- Allow for use of Monarch Hand Grips
- Cutout area (VLWGB05C) allows for use of thermoplastic and head support



VLWGB04C



VLWGB05C with Head Support Cutout

## Product Information

Ref #	Product Description
VLWGB04C	Wing Board Vac-Lok, 108 x 55 cm, 17 L, Nylon
VLWGB05C	Wing Board w/ Head Support Cutout, 108 x 55 cm, 14 L, Nylon



### Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

COPYRIGHT © 2020, CIVCO, AND POSIFIX ARE REGISTERED TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. MONARCH, VAC-LOK AND WING BOARDS ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. PN - WINGBOARD-VAC-LOK Rev. A

**CIVCO**<sup>®</sup>  
Radiotherapy