

Clearly Invisible

Universal Couchtop™ 1 Piece

Load Rating

- Entire System: 200 kg / 440 lbs

Deflection

- The Universal Couchtop 1 Piece meets the deflection standards and requirements of IEC 60976(2007) Section 15.2.2 on the Varian Exact® couch

Attenuation

- 6MV

Gantry	Transmission (10 x 10)	Transmission (5 x 5)
180°	0.971	0.969
200°/160°	0.968	N/A
210°	0.967	N/A
220°/140°	0.962	N/A

Surface Dose

	Surface Dose			
	6MV - 10 x 10	6MV - 20 x 20	18MV - 10 x 10	18MV - 20 x 20
Open Field	19.1%	30.5%	16.1%	29.6%
UCT 1 Piece	94.1%	97.2%	69.9%	79.4%