

## C-Qual M™ Breastboard and/or Monarch™ Overhead Arm Positioner

Patient:

Date:


Notes:

### Head Support Used


**1** Index Lok-Bar™ to table (ex: H1, F2)

*If using Monarch with Lok-Bar only, skip to Step 5*


**2** Set Tilting Board Elevation Angle (5° - 25°)



**3** Position Bottom Stop (B0 - B10)



**4** Set Monarch Overhead Arm Positioner (P0 - P77)



**5** Position Hand Grip (G1 - G12)

Hand Placement	Top	Side
----------------	-----	------

**6** Laser Alignment Location

**7** Thermoplastics Used

Head-Only	Yes	No
Breast	Yes	No