

C-Qual™ Breastboard with Monarch™ Overhead Arm Positioner

Patient:



Date:

Notes:


Head Support Used

1 Index Lok-Bar™ to table (ex: H1, F2)


2 Set Tilting Board Elevation Angle (5° - 25°)



3 Position Bottom Stop (B0-B10)




4 Set Monarch™ Overhead Arm Positioner (P0-P77)



5 Position Hand Grip (G1-G12)

Hand Placement	Top	Side

6 Laser Alignment Location



7 Thermoplastics Used

Head-Only	Yes	No
Breast	Yes	No