


Body Pro-Lok ONE™ ONEPlatform™

CIVCO RT Technical Data Sheet




Product Name	ONEPlatform™
Product Number	BPL003
Kits that include this product	BPL100 BPL110
General Description	An integral part of the Body Pro-Lok ONE System, the Platform features the same profile and indexing options as CIVCO's Universal Couchtop, providing seamless integration of positioning devices via three-pin Lok-Bar indexing.
Indications for use	The device is intended to aid in support and positioning of patients during MR, radiological and other procedures.
Contraindications	None
MRI Compatible	Device is MR safe.
Integration	Connects via two or three-pin Lok-Bar.
FDA Clearance	Yes
CE Marked (EU DOC)	Yes
Other Country Registrations	Available upon request
Country of Manufacture	USA

Technical Details

Primary Materials	Kevlar, Foam Core, Phenolic Inserts		
Weight (approx)	7.5 kg / 16.5 lbs		
Dimensions (mm)	1829 L x 535 W x 35 H (min.)		
Dosimetric Properties	Location	Field Size	Attenuation 6 (MV)
	Center of Platform	10 cm x 10 cm	2.1%

Technical Drawing



Packaging

Unit Pack	1 Piece
Units / Box	1 Piece
Shelf Life	Not Applicable
Warranty	See CIVCO RT website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	September 2014

Accessories

System Customization	ONEBridge™: BPL001 Standard Respiratory Plate: MTSBRT003 Wide Respiratory Plate: BPL008 Clam-Lok™ Cushion: MTSBRT202 Body Pro-Lok ONE™ Respiratory Belt: BPL026 Patient Handle, right: MTIL662501 Patient Handle, left: MTIL662502 Shoulder Restraint Bridge 2: MTSBRT038 Forehead Restraint: MTSBRT018 T-Shape Vac-Lok™ Cushion: MTVL0017 Monarch Overhead Arm Positioner: OPA110 Extended Wing Board w/ T-Grip: MTWB08T Extended Wing Board w/U-Grip: MTWB08U 3-Pin Lok-Bar™, Prodigy 2, 53 cm: MTIL6054 Slide Guide Lok-Bar: MTVL00853 Transverse Lok-Bar, Prodigy 2, 53 cm: MTIL6051 Kneefix™ 3 Cushion: 301052 Feetfix™ 3 Cushion: 301071
----------------------	---

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com



COPYRIGHT © 2022. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. BODY PRO-LOK ONE, CIVCO RADIOTHERAPY, CIVCO RT, ONEBRIDGE, CLAM-LOK, UNIVERSAL COUCHTOP, ONEPLATFORM, VAC-LOK, LOK-BAR, KNEEFIX, FEETFIX, WING BOARD, AND MONARCH ARE ALL TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
2020C1518 REV. E