

Body Pro-Lok™

CIVCO Technical Data Sheet



Product Name	Body Pro-Lok™ System	
Product Number	Sold in kits	
Kits Options	Full Platform Package Carbon Fiber: MTSBRT001	Kevlar: MTSBRT400K
	Lite Platform Package Carbon Fiber: MTSBRT001L	Kevlar: MTSBRT410K
	Ultra-Lite Platform Package Carbon Fiber: MTSBRT320	Kevlar: MTSBRT420K
	Comprehensive Body Pro-Lok UCT Package: MTSBRT310	
	Respiratory Restriction Body Pro-Lok UCT Package: MTSBRT300	
General Description	The modular design of the Body Pro-Lok System and accessories aids in providing comfortable, total body immobilization for lengthy hypo-fraction treatments.	
Indications for use	SBRT Patient Immobilization	
Contraindications	None	
MRI Compatible	MR Safe components are available. See product labeling or contact CIVCO for additional details.	
Integration	<ul style="list-style-type: none"> • Connects to most couchtops via standard Lok-Bar™ • Bridges connect to Universal Couchtop with rail adapters • All other accessories connect via 2-pin indexing 	
FDA Clearance	Yes	
CE Marked (EU DOC)	Yes	
Other Country Registrations	Available upon request	


Technical Details

Primary Materials Acetal (Black), Carbon Fiber (Platform)

Weight (approx) 19.5 kg / 43 lbs

Dimensions (mm) 2108 L x 637 W x 440 H (min.)

Technical Drawing



Country of Manufacture USA

Packaging

Unit Pack	1 Piece
Units / Box	1 Piece
Shelf Life	Not Applicable
Warranty	See CIVCO website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	July 2011

Accessories

System Customization	Respiratory Plate: MTSBRT003 Respiratory Belt: MTSBRT004 Type 1B Bridge: MTSBRT039 Type 2 Bridge: MTSBRT007 Type 3 Bridge: MTSBRT031 Patient Handle: MTSBRT005 Patient Handle, mirror: MTSBRT0051 Shoulder Restraint Bridge 2: MTSBRT038 Forehead Restraint: MTSBRT018 Clam-Lok™ Cushion: MTSBRT202 T-Shape Vac-Lok™ Cushion: MTVL0017 T-Pin Lok-Bar™: MTSBRT009 Kneefix™ 3 Cushion: 301052 Feetfix™ 3 Cushion: 301071 Transfer System: MTSBRT011 Type 2 & 3 Bridge Adapter, Qty 2: MTIL6615 Type 1B & Shoulder Restraining Bridge 2 Adapter, Qty 2: MTIL6616 Respiratory Belt Adapter, Qty 2: MTIL6620 Storage Cart: MTSBRT203 Wall Storage Unit: MTSBRT224
-----------------------------	---

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

