Bolus CIVCO RT Technical Data Sheet



Bolus with Skin, Skinless Bolus		
Skinless: 32052, 32102, 32152, 32202		
With Skin: 31051, 31101, 31151 31152, 31102, 31052, 31202, 31053050, 31103050, 31054040, 31104040, 31033030		
Bolus is made from a solid, homogeneous, tissue-equivalent gel. It will not sag, deform or run. Bolus is available encased in a tough layer of thin plastic "skin" or "skinless", which is sticky for greater conformity to steep body contours.		
Designed as bolus material for use in radiation therapy, diagnostic x-ray, and ultrasound imaging procedures.		
No		
Yes		
Yes		
Available upon request		
USA		



Technical Details

Primary Materials	Solid, homogenous, tissue-equivalent gel. Skin: Polyurethane/Synthetic Rubber			
Weight (approx) + Dimensions	Part Number		Weight (gm)	Dimensions
	31033030	With Skin	278.1	30 x 30 x 0.3 cm
	31053050	With Skin	772.5	30 x 50 x 0.5 cm
	31103050	With Skin	1545	30 x 50 x 1.0 cm
	31054040	With Skin	824	40 x 40 x 0.5 cm
	31104040	With Skin	1648	40 x 40 x 1.0 cm
	32052	Skinless	748.4	30 x 30 x 0.5 cm
	31052	With Skin	463.5	30 x 30 x 0.5 cm
	32102	Skinless	918	30 x 30 x 1.0 cm
	31102	With Skin	927	30 x 30 x 1.0 cm
	32152	Skinless	2177.2	30 x 30 x 1.5 cm
	31152	With Skin	1390.5	30 x 30 x 1.5 cm
	32202	Skinless	2812.3	30 x 30 x 2.0 cm
	31202	With Skin	1854	30 x 30 x 2.0 cm
	31051	With Skin	115.9	15 x 15 x 0.5 cm
	31101	With Skin	231.8	15 x 15 x 1.0 cm
	31151	With Skin	353	15 x 15 x 1.5 cm

Packaging

Unit Pack	1 piece
Units / Box	1 piece
Shelf Life	n/a
Warranty	See CIVCO RT website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	May 2006

Global Sales Office Orange City, Iowa, USA 800.842.8688 | +1 712.737.8688

info@CivcoRT.com www.CivcoRT.com



COPYRIGHT © 2022. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2017P1300 REV. B