

Chabner XRT® Radiation Bra

CIVCO Technical Data Sheet



Product Name	Chabner XRT® Radiation Bra		
Product Number			
	Individual Size	Radiation Bra Qty 1	Radiation Bra Qty 5
	1	CRTBRA01	CRTBRA01PK5
	2	CRTBRA02	CRTBRA02PK5
	3	CRTBRA03	CRTBRA03PK5
	4	CRTBRA04	CRTBRA04PK5
	5	CRTBRA05	CRTBRA05PK5
	6	CRTBRA06	CRTBRA06PK5
	7	CRTBRA07	CRTBRA07PK5
	8	CRTBRA08	CRTBRA08PK5
	9	CRTBRA09	CRTBRA09PK5
Kits that include this product	CRTBRASPK1-5 CRTBRAKIT1-5	CRTBRASPK5-9 CRTBRAKIT5-9	CRTBRASPK5-9 CRTBRAKIT1-9
General Description	Provides optimal breast support during simulation, planning and radiation treatment		
Indications for use	This device is indicated to aid in supporting adolescent and adult patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning. The device is not intended for use with patients under 12 years of age.		
MRI Compatible	Yes		
FDA Clearance	Yes		
CE Marked (EU DOC)	Yes		
Other Country Registrations	Available upon request		
Country of Manufacture	USA		

Technical Details

Primary Materials


UBL Tricot, Mesh Warp Knit, Synthetic Elastomer, Nylon, Polyester

Dosimetric Properties

Material	6 MV photons	6 MeV electrons	Proton water equivalency R80
Tricot	0%	0%	–
TPU ●	0%	0%	0.2 mm
Elastic	0%	.1%	–
Velcro ●	0%	.7%	–
Mesh ●	0%	0%	0.1 mm

Performance Specifications are for reference only, values may vary based on setup and equipment.

Technical Drawing



Part Number	Dimensions: Total Bottom Band on the Half (Fully Closed)	Product Weight (Approx)
CRTBRA01	14.25 in (36.2 cm)	2.1 oz (60 grams)
CRTBRA02	14.25 in (36.2 cm)	2.3 oz (65 grams)
CRTBRA03	15.75 in (40 cm)	2.5 oz (70 grams)
CRTBRA04	15.75 in (40 cm)	2.6 oz (75 grams)
CRTBRA05	16.00 in (40.6 cm)	2.8 oz (80 grams)
CRTBRA06	17.5 in (44.5 cm)	3.0 oz (85 grams)
CRTBRA07	17.375 in (44.1 cm)	3.2 oz (90 grams)
CRTBRA08	18.875 in (47.9 cm)	3.4 oz (95 grams)
CRTBRA09	18.75 in (47.6 cm)	3.5 oz (100 grams)

Packaging

Unit Pack	Each
Units / Box	One device
Shelf Life	n/a
Warranty	See CIVCO website for Product Warranty Statement
Sterilization	Non-Sterile
Launch Date	October 2017 – USA March 2018 – Global

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

