## Combifix<sup>™</sup> 3

### CIVCO Technical Data Sheet









Product Name	Combifix 3	
Product Number	151004MR	Combifix 3 for lower extremity positioning. Includes Combifix 3 Baseplate, Kneefix™ 3 Cushion, Elevation Block 3, Feetfix™ 3 Cushion and Feetfix 3 Couch Support.
	151005MR	Combifix 3 for lower extremity positioning. Includes Combifix 3 Baseplate, Low Kneefix 2 Cushion, Elevation Block 3, Feetfix 3 Cushion and Feetfix 3 Couch Support.
	302071MR	Feetfix 3 Couch Support (cushion not included)
General Description	The Combifix 3 is a baseplate designed to combine the Kneefix 3 and Feetfix 3 cushions, providing comfortable positioning with enhanced immobilization when treating the pelvic region and other lower extremities.	
Indications for use	The device is intended for positioning the hip and lower extremities of a patient during radiotherapy and diagnostics and is used to aid in the support and positioning of patients during an MRI.	
MRI Compatible	Yes	
Integration	2 or 3-Pin indexable	
FDA Clearance	Yes	
CE Marked (EU DOC)	Yes	
Other Country Registrations	Available upon request	





Country	of Manufacture	
---------	----------------	--


USA


# Technical Details

Product Name	Combifix 3	
Product Number	151004MR	
Primary Materials	PVC (Baseplate), Polyethylene (Cushions)	
Weight (approx)	20.3 lbs (9.2 kg)	
Dimensions (mm)	915 L x 580 W x 273 H	












### **Technical Drawing**





#### **Technical Drawing**

