

Medical Precision B.V. Comfort Marker 2.0

CIVCO Technical Data Sheet



Product Name	Comfort Marker 2.0	
Product Number	MPCM2010INT	Medical Precision Comfort Marker 2.0 System, Qty 1
	MPCM4022INT	Medical Precision Patient Marking Set. Includes: Safety Needle (Qty 50), Ink Cup (Qty 50), Black RadSafe Pigment (Qty 50)
	MPCM5022INT	Medical Precision Patient Marking Set. Includes: Safety Needle (Qty 50), Ink Cup (Qty 50), White RadSafe Pigment (Qty 50)
General Description	The Comfort Marker 2.0 allows the user to place accurate reference points on patients enabling radiotherapy. It consists of a Pen module which drives the Safety needle. Different depth settings can be chosen to accommodate for different skin types. When the reference point has been placed, the Control Unit can be stored and charged in the Docking station.	
Indications for use	Placement of markings on human skin at a controlled depth to enable effective radiotherapy.	
MRI Compatible	Not intended for use in MRI (magnetic resonance imaging) room	
FDA Clearance	No	
CE Marked (EU DOC)	Yes	
Other Country Registrations	Available upon request	
Country of Manufacture	The Netherlands	

Technical Details


Primary Materials Lithium polymer battery, Electronics, Metal and Plastics

Weight (approx) 0.7 kg (1.7 lb)

Dimensions 24.9 cm x 18.5 cm x 5.4 cm (9.8" x 7.3" x 2.1")

Unit/Box One



Technical Drawing (MPCM2010INT Base Unit)



Technical Details

Primary Materials	Zylar co-polymer (Housing); Stainless Steel (Spring/Needle)
Weight (approx)	1.19 kg (2.62 lb)
Dimensions	4.4 cm x 1.4 cm (1.7" x 0.5")
Quantity	Qty 50 (2 boxes of 25 individual pouches per box)
Shelf Life	Sterile, 1 year

Technical Drawing (Safety Needle)



Technical Details

Primary Materials	Biotic Phocea Radsafe Pigment
--------------------------	-------------------------------

Weight (approx)	3 g (0.006 lb)
------------------------	----------------

Dimensions	4.6 cm x 1.4 cm x 0.5 cm (1.8" x 0.5" x 0.2")
-------------------	---

Quantity	Pack of 50
-----------------	------------

Shelf Life	Sterile, 5 years
-------------------	------------------

Technical Drawing (Black or White Ink)



Technical Details

Primary Materials	Sterile Cup
Weight (approx)	1 gram
Dimensions	1.6 cm x 1.3 cm (0.6" x 0.5")
Quantity	Pack of 50
Shelf Life	Sterile, 3 years

Technical Drawing (Ink Cup)



Warranty

Warranty	See CIVCO website for Product Warranty Statement
Launch Date	August 2020

Global Sales Office
Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

