


HipFix® Baseplate

CIVCO Technical Data Sheet



Product Name	HipFix® Hip & Pelvic Immobilization System
Product Number	MTPLVC04MR
General Description	The HipFix System provides immobilization to enhance fixed fiducial targeting. The versatile system accommodates prone and supine setups while utilizing a baseplate and sheet of thermoplastic over the patient's abdomen and pelvic region.
Indications for use	Device is intended for patient fixation. Devices are used to aid in the support and positioning of patients during an MRI.
MRI Compatible	Yes
Integration	2 or 3-Pin Lok-Bar™
FDA Clearance	Yes
CE Marked (EU DOC)	Yes
Other Country Registrations	Available upon request
Country of Manufacture	USA


Technical Details

Product Name	HipFix®
Product Number	MTPPLVC04MR
Primary Materials	Black PVC
Weight (approx)	10.15 lbs (4.6 kg)
Dimensions (mm)	787 L x 603 W x 49 H (max)

Dosimetric Properties

Location	Transmission	Attenuation 6 (MV)	Attenuation 18 (MV)
Solid Region ●	10 cm x 10 cm	4.69%	3.31%

Technical Drawing



Packaging

Unit Pack	Each
Units / Box	One
Shelf Life	n/a
Warranty	See CIVCO website for Product Warranty Statement
Sterilization	Non-Sterile

Accessories

System Customization

Thermoplastics

Perforation	Handle Type	Thickness	Part #
Solid	Reloadable	3.2 mm	MTHFX1822S
Solid	Reloadable	2.4 mm	MTHFX18222.4S
Solid	Disposable*	3.2 mm	MTPLHA0201
Perforated	Reloadable	3.2 mm	MTHFX1822
Perforated	Reloadable	2.4 mm	MTHFX18222.4
Perforated	Disposable*	3.2 mm	MTPLHA0202
IMRT Reinforced	Reloadable	3.2 mm	MTHFI1822
IMRT Reinforced	Reloadable	2.4 mm	MTHFI18222.4
IMRT Reinforced	Disposable*	3.2 mm	MTPLHAI0203

*Groin piece sold separately

All disposable thermoplastics are MRI compatible

Handles & Groin Piece

MTPLHA01	Reloadable Handles (x2) & Groin Piece
MTPLVCSUB8	Groin Piece for use with disposable thermoplastics

Vac-Lok™ and RediFoam™

VL71FC	Split-Leg Vac-Lok, Nylon w/ indexing
MTFA39B	RediFoam, Small Includes 33 x 40" (83.8 x 101.6 cm) RediFoam bag
MTFA40B	RediFoam, Large Includes 39 x 46" (99.1 x 116.8 cm) RediFoam bag
MTFA41B	RediFoam, X-Large Includes 39 x 56" (99.1 x 142.2 cm) RediFoam bag

REDIFOAM ONLY AVAILABLE IN THE U.S.

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

