



Slide Guide Three-Pin Lok-Bar™

TECHNICAL DATA SHEET

PRODUCT NAMES	Slide Guide Three-Pin Lok-Bar
PRODUCT NUMBERS	MTIL6058
GENERAL DESCRIPTION	The Slide Guide Three-Pin Lok-Bar slides longitudinally along the Universal Couchtop™ Prodigy™ 2 profile without being limited to any indexing points. It may also be used to prevent lateral movement for longer devices such as bellyboards, other breastboards and Vac-Lok™ Cushions, which contract during forming.
INDICATIONS FOR USE	Devices are used to aid in the support and positioning of patients during an MRI.
MR SAFE	Yes
OTHER COMPATIBLE MODALITIES	CT, Photon, Proton
INTEGRATION	Three-Pin indexable to Universal Couchtop or TotalRT™ Platform profile.
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request.
COUNTRY OF MANUFACTURE	USA

PACKAGING

UNIT PACK	Each
UNITS / BOX	One device
SHELF LIFE	n/a
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Not sterile

TECHNICAL DETAILS

COMPONENT	MTIL6058 Slide Guide Three-Pin Lok-Bar
PRIMARY MATERIALS	Polyoxymethylene, Garolite G10, Stainless steel screws
WEIGHT (APPROX)	0.19 kg (0.4 lb)
DIMENSIONS (MM)	157 L x 607 W x 26 H
TECHNICAL DRAWING	

