

# Chin Strap Thermoplastic

The Posicast® Chin Strap Thermoplastic attaches to the Monarch Overhead Arm Positioner™ and any CIVCO RT Posifix® baseplate, extension or overlay.

## Features & Benefits

- Positions the patient's chin away from the treatment beam if treating nodes in the neck region
- Immobilizes the head in a reproducible position
- Drape style thermoplastic that requires less stretching which results in minimal shrinkage during cooling

772021 Chin Strap Thermoplastic, Posicast-Lite, Disposable, 2.3 mm, Qty 10



Monarch Overhead Arm Positioner



Posifix Baseplate



## Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

COPYRIGHT © 2022. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. POSICAST AND POSIFIX ARE REGISTERED TRADEMARKS OF MEDTEC LLC. CIVCO RADIO THERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2022P1635 REV A

**CIVCO**  
Radiotherapy™