

Head & Neck Immobilization for Proton Therapy

ProForm™

The Universal Couchtop™ ProForm™ Head & Neck Extension features geometry that allows for direct access to effectively treat CNS patients. A homogeneous treatment area and soft gradient edges reduce beam degradation. ProForm accepts a dedicated thermoplastic mask designed to meet the unique requirements of proton therapy.

Universal Couchtop™
Long Extension

No Overlay Required

- Large, 121.5 cm junction-free and artifact-free treatment area with quick and easy attachment of extension
- DoseMatch™ Technology – Dosimetric consistency via CT/SIM, photon & proton treatment compatible models
- Prodigy™ 2 indexing with 7 cm spacing
- Homogeneous, low density and low attenuation



Contact your sales representative to learn more about CIVCO's ProtonSeries™ family of solutions for proton therapy.

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

CIVCO®
Radiotherapy