

Suremark® Skin Markers

Suremark

Suremark markers are useful in distinguishing between a nipple shadow and lesion.



MTSL10	1.0 mm, Qty 110/Box
MTSL15	1.5 mm, Qty 110/Box
MTSL20	2.0 mm, Qty 110/Box
MTSL25	2.5 mm, Qty 110/Box

Suremark Tabs

Suremark Tabs maintain the advantages of the Suremark markers with an added tabbed feature.



MTSL15T	1.5 mm, Qty 110/Box
MTSL20T	2.0 mm, Qty 110/Box

Suremark Relief Tabs

Suremark Relief Tabs have a unique adhesive-free center that won't stick to sensitive areas of the skin.



MTSL25RT	2.5 mm, Qty 110/Box
----------	---------------------


Suremark Wire

Suremark Wire is used for marking scar tissue, outlining palpable nodules, moles and entry points in CT scans.



MTSLW3	0.3 mm, 336 cm
MTSLW4	0.4 mm, 336 cm
MTSLW8	0.8 mm, 336 cm
MTSLW15	1.5 mm, 336 cm

CT Mark



CT Mark new blue formulation is specifically designed for CT and will not cause artifact.

MTCTW10	Wire, 1.0 mm, 330 cm
MTCTW20	Wire, 2.0 mm, 300 cm
MTCT23	2.3 mm, Qty 110/Box
MTCT40	4.0 mm, Qty 50/Box
MTCTCR3	Cross, 1.2 cm, Qty 50/Box

Powermark™

Powermark markers maintain the advantages of Suremark markers and feature larger lead ball sizes for visibility in radiation therapy.



MTSL30	3.0 mm, Qty 50/Box
MTSL40	4.0 mm, Qty 50/Box

Clearline Scar Markers

Clearline Scar Markers provide outlining of scars and other areas of concern.



MTSLCL20	2.0 mm, 320 cm
MTSLCL25	2.5 mm, 203 cm

Visionmark

Visionmark markers are specially formulated non-metallic markers that perform consistently and accurately in a variety of applications.



MTSLV25	2.5 mm, Qty 110/Box
MTSLV30	3.0 mm, Qty 50/Box
MTSLV40	4.0 mm, Qty 50/Box
MTSLV50	5.0 mm, Qty 50/Box