# Universal Couchtop™ MRI Overlay & Adjustable Coil Holders

Siemens MAGNETOM® Sola®, Aera®, Vida® and Skyra®





A new dimension in radiotherapy treatment planning, enabling accurate transfer from MR imaging to treatment



The Universal Couchtop Overlay and Adjustable Coil Holders are designed to enable precise and reproducible patient positioning. The MRI Overlay is a lightweight board, made of fiberglass with a foam core. The overlay features Prodigy™ 2 indexing every 7 cm to help accurately transfer patient setup from imaging to treatment. Universal Couchtop is available on many major MR, CT, Proton and Linac models.

#### Standard Overlay Package

Universal Couchtop MR Overlay available on Siemens MAGNETOM Sola, Aera, Vida, Skyra MTM3004

(188 x 53 x 2.2 cm) with Type-S and Posifix® mask attachment. Includes Prodigy Three-Pin Lok-Bar™

(53 cm), Head Coil Support and Body Coil Support.

Accessory

MTM300404 Body Coil Support, Qty 1



# Positioning & Immobilization Device Compatibility

Universal Couchtop™ MRI Overlay attaches to compatible Siemens MAGNETOM® Sola®, Aera®, Vida® and Skyra® tables allowing for use of CIVCO RT patient positioning and immobilization devices during MR imaging for radiotherapy.

### **Head & Neck Positioning**

Compatible with both Type-S™ and Posifix® thermoplastics, including Head Only and Head, Neck & Shoulder options.



### **SBRT Positioning**

Body Pro-Lok ONE™ bridge and belt positioning devices secure directly to the overlay, allowing customized SBRT setups.



### Positioning with Three-Pin Lok-Bar™

A Three-Pin Lok-Bar allows indexing of CIVCO RT MR compatible devices to the overlay (example, Kneefix™ and Feetfix™). The three-pins help prevent non-MR compatible devices (two-pin) from being used on an MR system.



# Adjustable Coil Holders

## Positioning Flexibility When You Need It

Adjustable Coil Holders support and position the body and head coils, allowing the patient to be imaged in the same position as they will be in treatment.

#### Head

Customize the coil position with infinitely adjustable, "set and lock" coil holders.



## **Body**

Body coil support allows the coils to be adjusted for angled placement.





Orange City, Iowa, USA 800.842.8688 | +1 712.737.8688 info@CivcoRT.com www.CivcoRT.com

