
## TotalRT<sup>™</sup> All in One Platform TECHNICAL DATA SHEET



PRODUCT NAME	TotalRT All in One Platform
PRODUCT NUMBER	TRT100
GENERAL DESCRIPTION	The TotalRT All in One Platform features Prodigy <sup>™</sup> 2 indexing style and profile to quickly adapt many CQ Medical positioning and immobilization solutions to a patient's setup needs and treatment type. The platform is indexed on a couchtop creating an affordable all in one solution and allowing for seamless transition between patient treatments.
INDICATIONS FOR USE	The device is intended to aid in support and positioning of adult and pediatric patients undergoing radiation therapy including electron, photo, and proton treatments. The device is also used during image acquisition to support treatment planning.
MR SAFE	Yes
INTEGRATION	2 or 3-Pin indexable to couchtop
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

## **TECHNICAL DETAILS**


COMPONENT	TRT100 TotalRT All in One Platform		
PRIMARY MATERIALS	Glass Fiber, Foam Core		
WEIGHT (APPROX)	7.3 kg / 16 lbs		
MAX PATIENT WEIGHT	220 kg / 485 lbs		
DIMENSIONS (CM)	199 L x 61 W x 2.3 H (max)		
	LOCATION	ATTENUATION 6 MV PHOTONS 10 X 10 FIELD SIZE (CM <sup>2</sup> )	ATTENUATION 18 MV PHOTONS 10 X 10 FIELD SIZE (CM <sup>2</sup> )
DOSIMETRIC PROPERTIES	Middle of platform 🔴	2.23%	1.61%
PROPERTIES	Head area 🥚	1.76%	1.24%
	Performance Specifications are for reference only, values may vary based on setup and equipment.		
TECHNICAL DRAWING			
	แน่นารหนังออกไลและตั้งการที่เคาะกำหารกับและได้และค่างการสัต		
minanimainainainainainaina	านต่าวแก่กรรมโตกหนึ่งการที่สามมีมาแต่กระหลังการที่สามหรื		4





info@CQmedical.com | CQmedical.com

©2024 CQ MEDICAL, ALL RIGHTS RESERVED. CQ MEDICAL, CIVCO RT ARE TRADEMARKS OF MEDTEC LLC ALTA AND QFIX ARE REGISTERED TRADEMARKS OF QFIX SYSTEMS LLC ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. <u>VIEW OUR COMPREHENSIVE STATEMENT</u> CQ2024T55 REV A





## PACKAGING

UNIT PACK	1 Piece
UNITS / BOX	1 Piece
SHELF LIFE	Not Applicable
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
RELEASE DATE	October 2021

## ACCESSORES

SYSTEM	CQ Medical positioning and immobilization devices attach directly to platform, or
CUSTOMIZATION	via 2 or 3-Pin Lok-Bars.



info@CQmedical.com | CQmedical.com