

COMING SOON!

A Modular Approach to Breast & Thorax Positioning



Evolve Breastboard

This innovative and lightweight solution allows quick and simple indexing of many Wing Boards™ which aids in reproducibility without sacrificing patient comfort. Many hand grips, headrests and thermoplastic options easily integrate to maximize positioning options based on patient or clinical need. The board is indexable to couchtops with or without an optional folding bottom board and features smooth adjustability between five elevation angles.


Integrated GlideTrack

Allows simple positioning of the Pro, Extended and Standard Wing Boards to accommodate patients of varied torso lengths



Smooth Adjustment

Support bar secures the board to five adjustable elevation angles




- CT, MR, photon and SGRT treatment compatible
- Indexing marks allow reproducible setup
- Side scales assist with alignment
- Compatible with breast thermoplastic
- Indexable to TotalRT Platform or other couchtops via two or three-pin Lok-Bar™



Pro Wing Board

- Contoured cutouts in shoulder region eliminate attenuation during posterior axillary boosts
- Pro Hand Grips easily move along track, allowing many positioning options for the patient's arms
- Four handle options can be used (Pro Handles, M Handle, T or U Grip)
- Two headrest options are indexable in eight positions
- Chin thermoplastic can attach to position patient's chin away from the beam



Optional Bottom Stop

Slant board is indexable to the couchtop with or without the optional bottom board



Evolve Pro Breastboard Full Package




Evolve Pro Breastboard

Includes Slant Board, Bottom Board, Bottom Stop, Pro Wing Board, Pro Handles and Gel Head Cup

MTM500

Breastboard Configurations (choose your Wing Board™ below)



Evolve Slant Board, Bottom Board and Bottom Stop


OAP240



Evolve Slant Board Only

OAP200

Wing Board Options




Pro Wing Board with Pro handles and Gel Head Rest, 59 cm

OAP220




Extended Wing Board with T-Grip handle, 57 cm
(compatible with Silverman or Timo headrest, sold separately)

MTWB08T



Extended Wing Board with U-Grip handle, 57 cm
(compatible with Silverman or Timo headrest, sold separately)

MTWB08U



Standard Wing Board with Delrin handles, 57 cm
(compatible with Silverman or Timo headrest, sold separately)

MTWB09

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com