

Confident Proton Therapy without Compromise

 **Universal Couchtop™**
Long Extension

 **ProForm™**

The Universal Couchtop™ provides confidence in your care with no rails or junctions in the treatment area to facilitate better access, delivery of dose and target coverage.

- Large, artifact free treatment area with quick and easy attachment of extension
- Extension features geometry that allows for direct access to effectively treat CNS patients
- Homogenous treatment area and soft gradient edges reduce beam degradation
- Prodigy™ 2 indexing with 7 cm spacing
- Thermoplastic mask attaches with a smaller profile frame and pin similar to CIVCO RT's Type-S™ system



Pediatric ProForm™ Extension

COMING SOON!



CIVCO
Radiotherapy™

ProForm™ Head & Neck Overlays

Standard Overlays

CT & Photon
Treatment Overlay



MRI Overlay
 MRSeries™

Pediatric Overlays

COMING SOON!

CT & Photon
Treatment Overlay



MRI Overlay
 MRSeries™

Thermoplastic Styles

FreedomView™

Head Only or Head, Neck & Shoulders
Standard or Pediatric sizes



IMRT Reinforced™ with Hair Openings

Head Only
Standard size



ClearVision™ 2

Head Only or Head, Neck & Shoulders
Standard or Pediatric sizes



Standard Perforated

Head Only
Standard or Pediatric sizes



IMRT Reinforced™

Head Only or Head, Neck & Shoulders
Standard or Pediatric sizes



Solid with Perforated Nose

Head Only
Standard size



info@CivcoRT.com | www.CivcoRT.com

Visit CivcoRT.com for more information!

CIVCO
Radiotherapy™