


Clearly Visible

Gold Fiducial Markers

Planning • Intrafractional Tumor Tracking

Pretreatment Imaging • Shifting

Implantable fiducial markers provide a highly effective method of ensuring accurate target localization for tumors/organs which move in respect to external anatomy.



	Gold Marker	Needle
Knurled Gold Markers (pre-loaded in sterile needle)	0.9 x 3.0mm	18ga x 12cm
	0.9 x 3.0mm	18ga x 20cm
	1.2 x 3.0mm	17ga x 20cm
	1.2 x 3.0mm	17ga x 30cm
	1.6 x 3.0mm	14ga x 20cm
	1.0 x 5.0mm	18ga x 20cm
	0.8 x 10.0mm	18ga x 20cm
	1.2 x 10.0mm	17ga x 20cm
Coupled Markers (pre-loaded in sterile needle)	0.8 x 3.0mm	18ga x 20cm
	1.0 x 3.0mm	17ga x 20cm
	1.0 x 3.0mm & 1.2 x 3.0mm	17ga x 20cm
Knurled Gold Markers (packaged loose)	0.8 x 3.0mm	
	0.9 x 3.0mm	
	1.2 x 3.0mm	
	1.0 x 5.0mm	

Special Features

- Knurled surface to promote stability
- Clearly visible in kV and MV
- Approved for all soft tissue sites
- Available as Coupled Markers and in Cyberknife sizes

Procedure Flexibility

- Available preloaded and ready to deploy
- Loose gold for alternative delivery methods
- Delivered ready for sterile procedures



WWW.CIVCO.COM/RO | 800.842.8688 | +1 712.737.8688



Radiation Oncology

COPYRIGHT © 2013. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2013A0893