

CIVCO Set-up Sheet

Body Pro-Lok ONE™ System with ONEPlatform™ or Universal Couchtop™


Patient:

Date:

Notes:

Important instructions, please read:

- When indexing to the Universal Couchtop, use the indexing marks along the side of the couchtop for indexing Lok-Bars™ and the indexing marks along the top of the couchtop for indexing Body Pro-Lok ONE devices.
- When transferring a setup from the Universal Couchtop™ or ONEPlatform™ to the Body Pro-Lok ONE Rails-Only System, the ONEBridge™ must be adjusted and set 50 mm lower than on the Universal Couchtop or ONEPlatform*.
- When transferring a setup from the Universal Couchtop or ONEPlatform™ to the iBEAM® evo Couchtop with the Body Pro-Lok ONE evo Adapter, reference Appendix A for available setup regions*.



*Separate set-up sheets for Body Pro-Lok ONE Rails-Only and evo Adapter are available.

If utilizing the ONEPlatform, start at #1


If indexing directly to the Universal Couchtop, skip to #4

1 Lok-Bars™ for indexing the ONEPlatform are indexed to table at:


2 Lateral Lok-Bar Position: F G H

3 Lateral Position (only if using a Transverse Lok-Bar)

A) Head end indexed at:



B) Foot end indexed at:




4 Respiratory Plate (check all that apply):


Standard

Wide


ComfortCare™ Cushion




- A) ONEBridge™ Side A:
- B) ONEBridge Side B:
- C) Docking Index:



- D) ONEBridge Lateral Position:




- E) Respiratory Plate Height:
- F) Respiratory Plate Lateral Position: D C E



- G) Set pressure between: and
(only if using ComfortCare Cushion)

5 Respiratory Belt (check if utilized)




- A) Docking Index:
- B) Left velcro strap indexed at: -
- C) Right velcro strap indexed at: -
- D) Set pressure between: and


6 Patient Hand Grips (check if utilized)

- Right hand grip indexed at:
- Left hand grip indexed at:

7 Clam-Lok™ Cushion (check if utilized)




- A) ONEBridge™ Side A:
- B) ONEBridge Side B:
- C) Docking Index:




- D) ONEBridge Lateral Position:


8 Forehead Brace (check if utilized)



- A) ONEBridge Side A:
- B) ONEBridge Side B:
- C) Docking Index:



- D) ONEBridge Lateral Position:



- E) Forehead Brace Height:
- F) Forehead Brace Lateral Position: D C E

9 Additional indexable accessory:
Indexed at:

10 Additional indexable accessory:
Indexed at:

11 Additional indexable accessory:
Indexed at:

12 Additional indexable accessory:
Indexed at:
