

CIVCO RT Set-up Sheet

C-Qual™ Breastboard with Arm & Wrist Supports

Patient:



Date:

Notes:


Head Support Used

1 Index Lok-Bar™ to table (ex: H1, F2)


2 Set Tilting Board Elevation Angle (5° - 25°)



3 Position Bottom Stop (B0-B10)



4 Head Support Position (H1-H5)




5 Upper Arm Support Rotation

Select one for each arm

Right (A-P)

Left (A-P)




6 Upper Arm Support Horizontal Position

Select one for each arm

Right (T-Z)

Left (T-Z)




7 Upper Arm Support Angle

Select one for each arm

Right (1-12)

Left (1-12)



CIVCO Set-up Sheet


C-Qual™ Breastboard

8 Lower Arm Support Height


Select one for each arm

Right (A1-A2)

Left (A1-A2)



9 Laser Alignment Location



10 Thermoplastics Used

Head-Only	Yes	No
Breast	Yes	No

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com