## CIVCO RT Set-up Sheet

### Combifix<sup>™</sup> 3


Patient: Date:

Notes:

Note: The Feetfix™ Couch Support can be utilized without the Combifix, if using the Feetfix Couch Support alone skip step 2.

1 Lok-Bars™ are indexed to table at (ex: H1, F2):

7 Feetfix Style Couch Support Position




Indexed to Combifix at (select one):

A B C D E F G

Н

3 Feetfix Style Couch Support Rotation



(select one):

5° 10° 15° 20°

4 Feetfix Cushion Style





(select one):

Feetfix 2 (red)

Feetfix 3 (champange)

5 Feetfix Bracket Style





(select one):

Feetfix 2 Bracket (used only with red Feetfix 2 cushion) Feetfix 3 Bracket (used only with champagne Feetfix 3 cushion)


Cont'd >



# CIVCO RT Set-up Sheet

## Combifix<sup>™</sup> 3

### 6 Kneefix™ Style & Number of Elevation Blocks



Kneefix Style (select one):


Standard (champagne or red)

Low (black)

Number of Elevation Blocks (select one):

0 1 2

#### 7 Shoulder Retractor Settings (if used)



Shoulder Retractor Used (select one):

Standard (red handle) Short (black handle)

Shoulder Retractor Set at (selet one):

A B C D E

Global Sales Office Orange City, Iowa, USA 800.842.8688 | +1 712.737.8688

info@CivcoRT.com www.CivcoRT.com

