


# CIVCO RT Set-up Sheet



## ProForm™ Overlay Patient Handles & Universal Couchtop™ Patient Handles

For use when transferring a setup from the ProForm Overlay to the Universal Couchtop ProForm Extension.




Patient ID:

Date:

Notes:

Overlay	Extension
E0	G2
E1	G3
E2	G4
E3	G5
E4	G6
E5	G7
E6	G8
E7	G9
E8	H0
E9	H1
F0	H2
F1	H3
F2	H4
F3	H5
F4	H6
F5	H7
F6	H8
F7	H9
F8	I0
F9	I1
G0	I2
G1	I3
G2	I4
G3	I5
G4	I6
G5	I7
G6	I8
G7	I9
G8	J0
G9	J1
H0	J2
H1	J3
H2	J4
H3	J5
H4	J6
H5	J7
H6	J8
H7	J9
H8	K0
H9	K1
I0	K2
I1	K3
I2	K4
I3	K5
I4	K6
I5	K7

### 1 Overlay and Extension Indexing:



#### Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

