


# MR Visible Fiducial Markers



## FusionCoil™ Markers

The FusionCoil marker images equally well in MR and CT. It consists of an open spiral gold coil with a special biocompatible alloy core that appears as a solid black dot in MRI. This is the first of its kind providing superior reference for image fusion as well as precise IGRT. FusionCoil Cartridge Kits are available to reduce the number of needle passes and minimize the chance of Edema or Pneumothorax.




FusionCoil in CT (left) and MR (right)



## PolyMark™ Markers

The PolyMark Marker is the first biocompatible polymer-based marker. It minimizes artifact, is visible using kV based IGRT solutions, MRI and is optimized for Proton Therapy.



PolyMark in CT (left) and MR (right)



### Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

COPYRIGHT © 2020. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. FUSIONCOIL, MRSERIES AND POLYMARK ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2016P1101 REV. E

